

Result view
CLART FluAVir

AT code:
 40515
 Rev. 15.1

Sample reference: 5
 Array ID: 000000040515 (E1)
 Analysis type: tmb end point detection
 Date and time: Tue Aug 18 18:16:13 2009

Virus

Virus	Result	Controls
New Influenza A H1N1	POSITIVE	Passed

Virus	Result	Controls
Influenza A	Negative	Passed
Influenza A H1N1	Negative	Passed
Influenza A H3N2	Negative	Passed

- This diagnostic kit detects and genotypes the following Influenza subtypes in human samples:

INFLUENZA A (SUBTYPES)
Seasonal H3N2 *
Seasonal H1N1 *
New Influenza A (H1N1/09)
Generic Influenza A *

* In certain cases, the test may not be able to discriminate between seasonal subtypes (H1N1 and H3N2) of Influenza A. These genotype results will be reported as 'non-typeable'.

** *In silico* analysis has determined that this diagnostic kit may detect as Generic Influenza A the following subtypes: seasonal H1N1, seasonal H3N2, H1N2, H7N7 and H9N2.