

SAFETY DATA SHEET

According to Regulation 1907/2006/CE

Doc. N° 402156 rev.2
Date of issue: 02/03/2015
Date of revision: 09/10/2018

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product identifier

Product name : **TBX Agar (Tryptone Bile X-gluc Agar)**
Product Number : **402156**

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, culture medium for microbiology

1.3 Details of the supplier of the safety data sheet

Company : Biolife Italiana S.r.l.
Viale Monza 272, 20128 Milano Italia
Tel : 0039 02 252091
Fax: 0039 02 2576428
E-mail: mktg@biolifeitaliana.it

1.4 Emergency telephone number

Emergency Phone : 0039 02-6610-1029 (Centro Antiveleni Niguarda Ca'
Granda- Milano)

2. HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture

Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008 [CLP]

Pictogram none
Signal word none
Hazard statement(s) none
Precautionary statement(s) none

Supplemental Hazard Statements none

Caution - this mixture contains a substance not yet fully tested.

2.3 Other hazards none

3. COMPOSITION/INFORMATION ON INGREDIENTS

3.2 Mixtures

Component	Classification	Concentration
Sodium 3-α,12-α-dihydroxy-5-β-cholan-24-oate – Sodio Desossicolato		
CAS-No 302-95-4	Acute Tox. 4; STOT SE 3;	< 10 %
CE-No 206-132-7	H302, H335	

For the full text of the H-Statements mentioned in this Section, see Section 16

4. FIRST AID MEASURES

4.1 Description of first aid measures

Biolife Italiana S.r.l., Viale Monza 272, 20128 Milan, Italy. Tel. n° ++39 02 25209.1, Fax n° ++39 02 2576428
E-mail: mktg@biolifeitaliana.it; Web Site: www.biolifeitaliana.it

If inhaled	If breathed in, move person into fresh air. If not breathing, give artificial respiration.
In case of skin contact	Wash off with soap and plenty of water.
In case of eye contact	Flush eyes with water as a precaution.
If swallowed	Never give anything by mouth to an unconscious person. Rinse mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

4.3 Indication of any immediate medical attention and special treatment needed

no data available

5. FIREFIGHTING MEASURES**5.1 Extinguishing media****Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

5.2 Special hazards arising from the substance or mixture

Carbon oxides, Sulfur oxides, Sodium oxides

5.3 Advice for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

5.4 Further information

no data available

6. ACCIDENTAL RELEASE MEASURES**6.1 Personal precautions, protective equipment and emergency procedures**

Avoid dust formation. Avoid breathing vapors, mist or gas. Ensure adequate ventilation.

6.2 Environmental precautions

Do not let product enter drains.

6.3 Methods and materials for containment and cleaning up

Sweep up and shovel. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For disposal see section 13.

7. HANDLING AND STORAGE**7.1 Precautions for safe handling**

Provide appropriate exhaust ventilation at places where dust is formed.

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place. Store in cool place. Hygroscopic. Avoid moisture.

7.3 Specific end uses

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Biolife Italiana S.r.l., Viale Monza 272, 20128 Milan, Italy. Tel. n° ++39 02 25209.1, Fax n° ++39 02 2576428
E-mail: mktg@biolifeitaliana.it; Web Site: www.biolifeitaliana.it

8.1 Control parameters

Components with workplace control parameters

It doesn't contain substances with occupational exposure limit value.

8.2 Exposure controls

Appropriate engineering controls

Handle in accordance with good industrial hygiene and safety practice.

Personal protective equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.
The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Body Protection

The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

a) Appearance	beige powder
b) Odour	no data available
c) Odour Threshold	no data available
d) pH	7.0 – 7.4 at 25 °C at use concentration
e) Melting point/freezing point	
Melting point/range:	no data available
f) Initial boiling point and boiling range	no data available
g) Flash point	no data available
h) Evaporation rate	no data available
i) Flammability (solid, gas)	no data available
j) Upper/lower flammability or explosive limits	no data available
k) Vapour pressure	no data available
l) Vapour density	no data available
m) Relative density	no data available
n) Water solubility	no data available
o) Partition coefficient: noctanol/water	no data available

p) Auto-ignition temperature	no data available
q) Decomposition temperature	no data available
r) Viscosity	no data available
s) Explosive properties	no data available
t) Oxidizing properties	no data available

9.2 Other safety information no data available

10. STABILITY AND REACTIVITY

10.1 Reactivity	no data available
10.2 Chemical stability	no data available
10.3 Possibility of hazardous reactions	no data available
10.4 Conditions to avoid	no data available
10.5 Incompatible materials	Strong oxidizing agents
10.6 Hazardous decomposition products	no data available

11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects

Acute toxicity	no data available
Skin corrosion/irritation	no data available
Serious eye damage/eye irritation	no data available
Respiratory or skin sensitisation	no data available
Germ cell mutagenicity	no data available
Carcinogenicity	IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
Reproductive toxicity	no data available
Specific target organ toxicity - single exposure	no data available
Specific target organ toxicity - repeated exposure	no data available
Aspiration hazard	no data available
Potential health effects	
Inhalation	May be harmful if inhaled. Causes respiratory tract irritation.
Ingestion	May be harmful if swallowed.
Skin	May be harmful if absorbed through skin. Causes skin irritation.

Eyes	Causes serious eye irritation.
Signs and Symptoms of Exposure	To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.
Additional Information	RTECS: no data available

12. ECOLOGICAL INFORMATION

12.1 Toxicity	no data available
12.2 Persistence and degradability	no data available
12.3 Bioaccumulative potential	no data available
12.4 Mobility in soil	no data available
12.5 Results of PBT and vPvB assessment	no data available
12.6 Other adverse effects	no data available

13. DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods	
Product	Offer surplus and non-recyclable solutions to a licensed disposal company.
Contaminated packaging	Dispose of as unused product.

14. TRANSPORT INFORMATION

14.1 UN number	ADR/RID: - IMDG: - IATA: -
14.2 UN proper shipping name	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.3 Transport hazard class(es)	ADR/RID: - IMDG: - IATA: -
14.4 Packaging group	ADR/RID: - IMDG: - IATA: -
14.5 Environmental hazards	ADR/RID: no IMDG Marine pollutant: no IATA: no
14.6 Special precautions for user	no data available

15. REGULATORY INFORMATION

This safety datasheet complies with the requirements of Regulation (EC) No. 1907/2006.

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture	no data available
15.2 Chemical Safety Assessment	no data available

16. OTHER INFORMATION

Text of H-code(s) and R-phrases mentioned in Section 3

Acute Tox.	Acute toxicity
H302	Harmful if swallowed.

H335
STOT SE

May cause respiratory irritation.
Specific target organ toxicity - single exposure

Further information
Restrictions for use

No data available

Training advice

No data available

References

No data available

The information in this MSDS was obtained from current and reliable sources. Anyone using this information is solely responsible for the accuracy and applicability of this information to a particular use or situation. Biolife Italiana S.r.l. does not in any way warrant or imply the applicability, viability or use of this information to any person or for use in any situation.